

















Wet-Gel Electrodes

Feta Med's wet gel electrodes utilize a non-irritating Potassium Chloride gel to provide fast signal response with low impedance. The LL, LRM, and LLE Series electrodes employ our patented lead-locking feature to secure the lead wires and prevent motion artifacts associated with movement of the snap on the pin; they are ideal for ambulatory procedures. All of our electrodes are Latex-free.

Name/Specs	Quantity	Product Picture
<ul style="list-style-type: none"> • LL Electrode Series • Patented Lead-Locking System • Wet Gel • 1/16" thick foam, 1-3/4" X 2-3/8" 	<p>1/pk, 50/bx, 500cs 3/pk, 45/bx, 450cs 4/pk, 600/cs 5/pk, 75/bx, 750/cs 7/pk, 1050/cs</p>	
<ul style="list-style-type: none"> • LL911 Electrode • Patented Lead-Locking System • Wet Gel • 1/16" thick foam, 1-3/4" X 2-3/8" • Made for EMS Services 	<p>3/pk, 45/bx, 450/cs</p>	
<ul style="list-style-type: none"> • LRM Electrode Series • Patented Lead-Locking System • Wet Gel • 1/16" thick foam, 1-3/4" X 2-1/2" 	<p>3/pk, 45/bx, 450/cs 5/pk, 75/bx, 750/cs</p>	
<ul style="list-style-type: none"> • LLEBL Electrode Series • Patented Lead-Locking System • Wet Gel • 1/16" thick foam, 1-1/2" X 2-1/4" 	<p>3/pk, 600/cs 5/pk, 550/cs 7/pk, 700/cs 30/bulk pk, 900/cs 10/pk, 900/cs</p>	
<ul style="list-style-type: none"> • LLEBX Electrode Series • Wet Gel • 1/16" thick foam, 1-1/2" X 2" 	<p>3/pk, 600/cs 5/pk, 550/cs 30/bulk pk, 900/cs 10/pk, 900/cs</p>	

<ul style="list-style-type: none"> ● P6 Electrode ● Wet Gel ● 1/16" thick foam, 7/8" X 1-3/4" 	<p>3/pk, 108/bx, 1080/cs</p>	
<ul style="list-style-type: none"> ● F1010 Electrode ● Wet Gel ● 1/16" thick foam, 7/8" X 1-3/4" 	<p>10/pk, 100/bx, 1000/cs</p>	
<ul style="list-style-type: none"> ● 200-19, W-60 Electrode ● Wet Gel ● White cloth type circular shape 2-3/8" in diameter 	<p>30/bulk pk, 1200/cs</p>	
<ul style="list-style-type: none"> ● LL1028 Electrode Series ● Wet Gel ● 1/16" thick white foam ● 1-1/2" Teardrop / Tab 	<p>5/pk, 75/bx, 750/cs</p>	
<ul style="list-style-type: none"> ● LL1036 Electrode Series ● Wet Gel ● 1/16" thick white foam ● 1-7/8" Teardrop / Tab 	<p>30/bulk pk, 600/cs</p>	
<ul style="list-style-type: none"> ● LL1136 Electrode ● Wet Gel ● 1/16" thick white foam ● 1-1/8" X 2-1/8" Rnd. Rect. / Tab 	<p>50/bulk pk, 1000/cs</p>	

<ul style="list-style-type: none"> ● LL5110V Electrode ● Wet Gel ● 1/16" thick white foam ● 1-5/8" X 2-1/8" Rect. / Tab 	<p>10/pk, 100/bx, 1000/cs</p>	
<ul style="list-style-type: none"> ● LL5105V Electrode ● Wet Gel ● 1/16" thick white foam ● 1-5/8" X 2-1/8" Rect. / Tab 	<p>5/pk, 75/bx, 750/cs</p>	
<ul style="list-style-type: none"> ● LL5103V Electrode ● Wet Gel ● 1/16" thick white foam ● 1-5/8" X 2-1/8" Rect. / Tab 	<p>3/pk, 60/bx, 600/cs</p>	
<ul style="list-style-type: none"> ● 200-21, CT-50 Electrode ● Wet Gel ● Clear Tape Vinyl ● 2" Round 	<p>30/pouch, 1200/cs</p>	
<ul style="list-style-type: none"> ● 200-09, FS-50C Electrode ● Wet Gel ● White foam ● 2" Round ● Radiotranslucent Carbon snap 	<p>30/pouch, 1200/cs</p>	

pk=pack, bx=box, cs=case